

## ROTORS

MAX. VOLUME	MAX. SPEED	RCF VALUE	ORDER NUMBER	ANGLE
6 x 250 ml	8,000 rpm	10,016xg	Z366-06250	28°
4 x 250 ml	4,500 rpm	3,780xg	Z366-250	0 – 90°
16 x 50 ml conical	4,500 rpm	3,735xg	Z366-250-HC	0 – 90°
6 x 85 ml	11,000/13,000 rpm	15,555/21,726xg	Z326-0685	38°
30 x 15 ml	4,500 rpm	2,830xg	Z326-3015	35 °
4 x 85 ml	12,000 rpm	14,809xg	Z326-0485	30°
6 x 50 ml round	13,000 rpm	15,869xg	Z326-0650-HS	26°
20 x 10 ml	12,000 rpm	15,775xg	Z326-2010	30°
44 x 1.5/2.0 ml	15,000/16,000 rpm	21,379/24,325xg	Z326-4420	32°/55°
24 x 1.5/2.0 ml	15,000/16,000 rpm	21,379/24,325xg	Z326-2420	45°
24 x 1.5/2.0 ml*	15,000/16,000 rpm	21,379/24,325xg	Z326-2420H	45°
10 x 50 ml conical	9,000/10,500 rpm	11,771/16,022xg	Z326-1050	38°
2 x 3 Microplates	4,500 rpm	2,716xg	Z326-02MP	0 – 90°

\* hermetically sealed

## EXAMPLE ROTORS

Main rotor



Max. volume: 4 x 250 ml  
Order number: Z366-250



Max. volume: 16 x 50 ml conical  
Order number: Z366-250-HC



Max. volume: 2 x 3 Microplates  
Order number: Z326-02MP



Max. volume: 6 x 250 ml  
Order number: Z366-06250



Max. volume: 6 x 85 ml  
Order number: Z326-0685



Max. volume: 10 x 50 ml conical  
Order number: Z326-1050

The respective max. capacities of the rotors are shown in the tables. Naturally, adapters are available for smaller volumes. Detailed information can be found in our individual brochures on our homepage: [www.HermleUSA.com/download](http://www.HermleUSA.com/download)  
If you cannot find the required adapter, we will gladly make a custom adapter for you.