

## ROTORS

MAX. VOLUME	MAX. SPEED	RCF VALUE	ORDER NUMBER	ANGLE
44 x 1.5/2.0 ml	15,000 rpm	21,127 / 21,379 xg	Z216-4420	32°/55°
30 x 1.5/2.0 ml	13,500 rpm	19,151 xg	Z216-3020H	45°
12 x 5 ml	14,500 rpm	19,978 xg	Z216-1205	45°
24 x 1.5/2.0 ml*	15,000 rpm	21,379 xg	Z216-2420H	45°
24 x 1.5/2.0 ml	15,000 rpm	21,379 xg	Z216-2420	45°
64 x 0.5 ml	13,500 rpm	14,261 / 16,298 xg	Z216-6405	40°
4 x 8-PCR strips	15,000 rpm	15,343 xg	Z216-3202	45°
24 Capillaries	12,000 rpm	14,970 xg	Z216-24HEM	90°
24 x 1.5/2.0 ml + 4 x 8-PCR	15,000 rpm	21,631 xg	Z216-2420S	45°
24 x 2.0 ml Spin Column Kits	15,000 rpm	21,379 xg	Z216-24SC	45°

\* hermetically sealed

## EXAMPLE ROTORS

Main rotor



Max. volume: 24 x 1.5/2.0 ml

Order number: Z216-2420



Max. volume: 12 x 5 ml

Order number: Z216-1205



Max. volume: 44 x 1.5/2.0 ml

Order number: Z216-4420



Max. volume: 24 x 1.5/2.0 ml

Order number: Z216-2420H

Hermetically sealed



Max. volume: 4 x 8-PCR strips

Order number: Z216-3202



Max. volume: 24 Hematocrit capillaries

Order number: Z216-24HEM *only Z216-M*

### Combi-Rotor



Max. volume: 24 x 1.5/2.0 ml, 4 x 8-PCR

Order number: Z216-2420S



Max. volume: 24 x 2.0 ml Spin Column Kits

Order number: Z216-24SC